




No Compromise

ACROSS NC

PTCA RX Catheter


- Real Non-Compliant catheter
- High RBP
- Fully Hydrophilic Coating
- Excellent Re-Wrap
- Controlled equal growth



Double marker

Low entry
profile 0.016"

High RBP up to 21 Atm



Ordering Information / Technical Information

Length (mm)	Diameter (mm)							
	2.00	2.50	2.75	3.00	3.25	3.50	4.00	4.50
5*	ACRO2005	ACRO2505	ACRO2705	ACRO3005	ACRO3205	ACRO3505	ACRO4005	ACRO4505
10	ACRO2010	ACRO2510	ACRO2710	ACRO3010	ACRO3210	ACRO3510	ACRO4010	ACRO4510
15	ACRO2015	ACRO2515	ACRO2715	ACRO3015	ACRO3215	ACRO3515	ACRO4015	ACRO4515
20	ACRO2020	ACRO2520	ACRO2720	ACRO3020	ACRO3220	ACRO3520	ACRO4020	ACRO4520

Order number description: ACRS DDLL (D=Diameter, L= length)

* ACROSS NC with 5 mm length has only 1 marker.

Control your Size



The balloon material is a NON - compliant polyamide for controlled dimensions

Pressure BAR	Atm 4	Atm 5	Atm 6	Atm 7	Atm 8	Atm 9	Atm 10	Atm 11	Atm 12	Atm 13	Atm 14	Atm 15	Atm 16	Atm 17	Atm 18	Atm 19	Atm 20	Atm 21
2.00 mm	1.74	1.78	1.82	1.85	1.9	1.93	1.96	1.98	2.00	2.03	2.06	2.09	2.11	2.14	2.16	2.18	2.21	2.24
2.50 mm	2.27	2.31	2.35	2.38	2.41	2.44	2.46	2.48	2.50	2.52	2.54	2.56	2.58	2.6	2.62	2.64	2.66	2.69
2.75 mm	2.46	2.48	2.51	2.55	2.6	2.65	2.68	2.72	2.75	2.77	2.79	2.82	2.84	2.86	2.89	2.91	2.93	2.96
3.00 mm	2.75	2.80	2.84	2.87	2.9	2.93	2.96	2.98	3.00	3.02	3.04	3.06	3.08	3.11	3.13	3.15	3.18	3.20
3.25 mm	2.83	2.86	2.91	2.96	3.01	3.08	3.13	3.18	3.25	3.27	3.3	3.33	3.37	3.41	3.44			3.48
3.50 mm	3.13	3.16	3.23	3.30	3.35	3.4	3.44	3.47	3.50	3.52	3.55	3.58	3.60	3.63	3.66			3.68
4.00 mm	3.57	3.65	3.73	3.80	3.86	3.91	3.95	3.98	4.00	4.03	4.06	4.08	4.11	4.14	4.17			4.20
4.50 mm	4.07	4.15	4.23	4.30	4.36	4.41	4.45	4.48	4.50	4.53	4.56	4.58	4.61	4.64	4.67			4.70

No Compromise

